



Strategies for Cleaning Process Validation in the Pharmaceutical Environment

Enabling automated cleaning of complex machine parts for an efficient validation

Melich Dietrich Seefeldt & Dr. Olaf Neuschaefer-Rube
20th April 2021

Outline

„Enabling automated cleaning of complex machine parts for an efficient validation“

- | | | |
|----------|---|----|
| 1 | Regulative definitions & automated cleaning | 6 |
| 2 | Complex parts' cleaning & validation | 9 |
| 3 | Customer trials & cases | 11 |

A man in a dark shirt and pants stands on a metal staircase inside a large, curved industrial structure, possibly a tunnel or a large pipe. He is looking up and reaching towards the ceiling. The structure is made of metal and has many pipes and cables running along its length. A red banner is overlaid across the middle of the image.

This is us – Fedegari

Company overview



- Fedegari founded in Pavia, Italy in **1953**
- Global provider of machines & integrated solutions for the **Contamination Control** in the **Bio-/ Pharmaceutical industry**
- Globally acting, family-owned group with 7 international subsidiaries and **550** employees world-wide



Fedegari portfolio

Turn-key solutions & required services for the techno-economic production of innovative sterile drugs





1. Regulative definitions & automated cleaning



Regulative definitions & automatic cleaning



Global regulatory bodies define cleaning validation similarly & state concerns about manual cleaning to deliver insufficient scientific data

Source

Definition validation

Automatic vs. manual application



FDA Guide to Inspections of Validation of Cleaning Processes (1993)

*„In the end, the test of any **validation** process is whether **scientific data** shows that the system **consistently does as expected** and produces a result that consistently **meets predetermined specifications**“*

„How variable are manual cleaning processes from batch to batch and product to product?“



Annex 15: Qualification and Validation (2015)

„Cleaning validation is documented evidence that an approved cleaning procedure will reproducibly remove the previous product or cleaning agents used in the equipment below the scientifically set maximum allowable carryover level.“

„10.15. Where manual cleaning of equipment is performed, it is especially important that the effectiveness of the manual process should be confirmed at a justified frequency.“



PIC/S 006-3 – Validation Master Plan (2007)

„7.1.3 ...Cleaning Validation is documented evidence that an approved cleaning procedure will provide equipment which is suitable for processing of pharmaceutical ingredients (APIs)“

„How many times need a manual cleaning process be applied to ensure adequate cleaning of each piece of equipment?“



Complex parts' cleaning & validation

Within the life cycle concept of the 2011 FDA guidance the design of the key-equipment-to-be-cleaned is an important factor, to reach a validatable cleaning procedure, but...

Stage 3: Process Monitoring

- Monitoring of Critical Process Parameters as Part of APR and Other Monitoring Programs



Stage 1: Process Design

- Define the Knowledge Space
- Identify Critical Process parameters
- Determine Control Strategy

„The **functionality** and limitations of commercial **manufacturing equipment** should be **considered** in the **process design**, as well as predicted **contributions to variability** posed by different component lots, production operators, environmental conditions, and measurement systems in the production setting.“

“...The results of **DOE studies** can provide **justification** for establishing **ranges** of incoming component quality, **equipment parameters**, and in-process material quality attributes...”

Stage 2: Process Qualification

- Equipment / Utility / Facility Qualification
- Process Performance Qualification



2. Complex parts' cleaning & validation



Complex machine parts

...many complex machine parts often cannot be cleaned automated, compared to standard ones


Complex machine parts

Examples:

- Tableting devices,
- Multi-port containers
- Custom-made production devices
- ...

Influencing factors on complexity:

- Accepted functional design,
- Residue acceptance criteria,
- Product quality expectations
- Production efficiency.
- ...



- Specifically **custom-engineered devices** in automatic washers can guarantee **full coverage** of most **complex load surfaces**.
- They can **treat** defined peculiarities and **enable standardization** and **proven reproducibility of procedures** and therefore the required “scientific data” for process validation.




The background of the slide is a black and white photograph of a shower stall. The shower door is closed, and the interior of the stall is visible through the glass. The shower floor has a grid pattern, and there are some dark spots on the wall. A red horizontal bar is overlaid across the middle of the image, containing the text.

3. Customer trials & cases



Reaching validated cleaning process through design & quality

Flexible & module-engineered design of washing racks for dedicated loads



The washing racks can include:

- Nozzles
- Spray arms
- Rotating spray balls
- Water blades
- Orifices










Automated equipment solutions for customer cleaning validation

Collaborative activities to facilitate verifiable cleaning strategies: Overview

User support by cleaning equipment supplier

User validation requirements

Prevention of cross-contamination		Sampling position	Analysis of best sampling position, detergent evaluation, equipment & process design
Predetermined residual levels		Detection systems	<u>Specific:</u> liquid sampling to monitor drug residual absence <u>Unspecific:</u> Conductivity meter for final rinse monitoring & 2 nd used also for detergent conc. determination, TOC available for measuring total organic residuals
Routine monitoring		R&D trials, rack design	Material testing by in-house R&D trials, engineering
Complexity of parts		R&D trials	Analysis of parts complexity and cleaning testing by in-house R&D trials
Material compatibility		Material, R&D trial, in-house production	Best quality material (surface 'BA' finish, 316L), in-house manufacturing of all critical components
Equipment design		Process design	With clean steam used as cleaning agent with emollient effect
Cleaning procedure		Automation system	Safe & reliable controller system: GAMP5, validated process recording system.
Documentation			



Case 1: Reaching cleaning requirements using design & quality

Collaborative activities to facilitate verifiable cleaning strategies: R&D trials

Removal of pharmaceutical drug (cytotoxic)

Manufacturer requirements: parts visual clean and dry

Identification of critical parts & areas for automated cleaning trials using Mannitol as 'placebo' contamination (shows similar properties as drug)



Identification of the worst part for automated cleaning






Case 1: Reaching cleaning requirements using design & quality

Collaborative activities to facilitate verifiable cleaning strategies: R&D trials

Removal of pharmaceutical drug (cytotoxic)

Trial activities simplify cleaning process verification tasks for routine cleaning process





Case 2: Reaching cleaning requirements using design & quality

Collaborative activities to facilitate verifiable cleaning strategies: R&D trials

Gelatine removal from capsule manufacturing part

**Manufacturer requirements: Reducing microbial contamination (< 30 cfu/ml)
and parts visual clean**

Assembled
production part



'Placebo' contamination simulation, using "coloured" gelatine for R&D cleaning trials



Artificially 'contaminated' production part
ready to perform 1st cleaning trials






Case 2: Reaching cleaning requirements using design & quality

Collaborative activities to facilitate verifiable cleaning strategies: R&D trials

Gelatine removal from capsule manufacturing part

Cleaning trials leads to engineered solutions towards automated cleaning process


Identifying critical areas during initial cleaning trials



Testing design: simultaneous steam injection (100°C) and mechanical action by air bubbling



Test resulting in 'visible clean' parts and microbial reduction is possible through subsequent sterilisation






Washer-Steriliser: Validated 'cleaning' & 'sterilization' in one unit

Automated cleaning & subsequent sterilization can be performed as validated processes with just a single equipment



FOWS Washer-Sterilizer features:

- Automated cleaning as validated process
- Saturated steam sterilization as validated process

- 
- Joint (and single) validated and automated cleaning & sterilization processes
 - Saving of installation space
 - Increased efficiency & time saving
 - Cost savings (in investment & maintenance)

Ihr Full-Service Partner für Sterilisationslösungen

Wir sind für Sie da:

Dr. Olaf Neuschäfer-Rube



- Technischer Vertrieb Industriemaschinen
- Prozessdesign
- Projektentwicklung



Melich Dietrich Seefeldt



- Leitung Niederlassung
- Vertriebs-, Angebots & Vertragsmanagement
- Investitionsberatung



Igor Vidari



- After-Sales zuständig
- Installation & start-up
- Validation Techniker



+49 (0)69 9074 66-200

Labor@de.fedegari.com



Industrie@de.fedegari.com



Ersatzteile@de.fedegari.com





Watch us on
YouTube



Thank you

www.fedegari.com

